

Pre-fabricated Knee Extension Brace

How to measure for a Pre-fabricated Knee Brace:


1. Position leg in maximum extension.
2. A circumference measurement should be taken around the leg at each level.

Level A - is located 9.5 inches (25 cm) above joint center.

Level B - is located at the joint center.

Level C - is located 9.5 inches (25 cm) below joint center.

Measurement for Pre-fabricated Knee Extension Brace			
		Circumference	
		Inches	Centimeters
Small	A	12.75 - 16.5	32.4 - 41.9
	B	11.75 - 14.5	29.8 - 36.9
	C	8 - 11.75	20.3 - 29.8
Medium	A	18 - 21	45.7 - 53.3
	B	15 - 17	38.1 - 43.2
	C	9.5 - 12.5	24.1 - 31.8
Large	A	21 - 24	53.3 - 61
	B	17 - 20	43.2 - 50.8
	C	12.75 - 16.5	32.4 - 41.9



Pre-fabricated Elbow Extension Brace

How to measure for a Pre-fabricated Elbow Brace:


1. Facing the patient, position the elbow in maximum extension and the wrist in neutral pronation/supination. See diagram to check proper positioning for measurement.
2. A circumference measurement should be taken around the arm at each level.

Level A - is located 7 inches (17.8 cm) above joint center.

Level B - is located at the joint center.

Level C - is located 7 inches (17.8 cm) below joint center.

Measurement for Pre-fabricated Elbow Extension Brace					
		Inches		Centimeters	
		circumference		circumference	
Small	A	8 - 11.75		A	20.3 - 29.8
	B	8 - 11.75		B	20.3 - 29.8
	C	5.75 - 8		C	14.6 - 20.3
Medium	A	10.25 - 14.25		A	26 - 36.2
	B	10.25 - 14.25		B	26 - 36.2
	C	6.75 - 10.5		C	17.1 - 26.7
Large	A	12.75 - 16.5		A	32.4 - 41.9
	B	12.75 - 16.5		B	32.4 - 41.9
	C	8 - 11.75		C	20.3 - 29.8




Pre-fabricated Wrist Extension Brace

How to measure for a Pre-fabricated Wrist Brace

1. See diagram to check proper positioning for measurements.

Level A - Circumference at forearm band/strap location.

Level B - M-L of hand across MP joints.



SMALL		Circumference	Diameter (M-L)
	LEVEL A	6 1/2" - 9"	
	LEVEL B		2 1/2" - 3 1/4"

MEDIUM	LEVEL A	9" - 11"	
	LEVEL B		3 1/4" - 4 1/4"

Hip Abduction Orthosis

How to measure for a Pre-fabricated Hip Abduction Orthosis:

1. Initial fitter is required to be proficient at bending metal. (BENDING IRONS REQUIRED, NOT INCLUDED)
2. Measure thigh circumference at levels "A" and "B" to determine size for thigh assembly.
3. Measure waist circumference. This product will fit a patient with a waist circumference up to 65 inches.

Thigh Assembly Sizing Chart		
Size	Prozimal (A)	Distal (B)
Small	15" - 22"	13" - 19"
Medium	18" - 25"	15" - 21"
Large	21" - 27"	17" - 23"

