

PCC #: _____

BILL TO: _____

ADDRESS: _____

SHIP TO: ☐ SAME AS BILLING _____

ADDRESS: _____

SHIPPING: ☐ GROUND (FXGD) ☐ STANDARD 2 DAY (FX2D)
OVERNIGHT: ☐ PRIORITY (FX1D) ☐ 1st OVERNIGHT (FX1A)
☐ OTHER: _____

CLINICIAN: _____

CELL #: _____

PATIENT ID: _____

HEIGHT: _____ **WEIGHT:** _____ **AGE:** _____

☐ MALE ☐ FEMALE
☐ LEFT ☐ RIGHT ☐ BILATERAL

NG ENCOUNTER #: _____

MEASUREMENT DATE: _____

IN-OFFICE REQUEST DATE & TIME: _____

☐ PROJECT NEWTON (Credit applied on prior auth. denial, details on One Hanger)

HiPro is not available digitally. Please ship cast directly to HFN Houston 3620 Willowbend Blvd, Ste 1118-A | Houston, TX 77054 **HFN:** ☐ HOUSTON ☐ OTHER _____

Check Brace Fitting & Casting Quick Guide

ACTIVITY LEVEL (Check one)

- ☐ Limited ambulator: sit to stand and transfer
- ☐ Household ambulator: level surfaces with walking aids
- ☐ Limited community ambulator: level surfaces with walking aids
- ☐ Active community ambulator: mild inclines and declines with or without walking aids
- ☐ Independent ambulator: varied cadence, uneven surfaces and no walking aids
- ☐ Active ambulator: walking, running, some athletic activity

OPTIONS

- ☐ Inner Boot ☐ Anterior Shell

CAST EVALUATION

- Side** **Heel Height**
- ☐ Left ☐ Right _____" ☐ Correct ☐ Leave as Cast
- Rotation**
- ☐ OK ☐ Int ☐ Ext _____° ☐ Correct ☐ Leave as Cast
- Fore Foot**
- ☐ OK ☐ Inv ☐ Evr _____" ☐ Correct ☐ Leave as Cast
- Rocker**
- ☐ OK ☐ Flat ☐ Contoured ☐ Correct ☐ Leave as Cast

GUIDANCE



- Any brace with a flat toe plate will require shoe work to function
- Any brace without a contoured foot plate/toe ramp is prone to fail

REQUIREMENTS

- Shoe insert or tracing
- Check orthosis fitting prior to final fabrication

Check List

- ☐ Patients shoe shipped with cast (preferred)
- ☐ Tracing of shoe insole provided with order form
- ☐ No reference provided (forefoot segment will be made large and will require adjustment by clinician)
- ☐ Use fiberglass for AFO section ☐ Use plaster for anterior shell on check orthosis



LAB USE ONLY

Spring Category (Strut length may change due to final design)

- ☐ .5 ☐ .75 ☐ 1 ☐ 2 ☐ 3 ☐ 4 ☐ 5 ☐ 6 ☐ 7

Spring Length

- ☐ 200mm ☐ 250mm ☐ 300mm

TURNAROUND TIMES

To review current projected turnaround times for fabrication sites visit the [Daily HFN Capacity Webpage](#)

PCC #: _____

BILL TO: _____

ADDRESS: _____

SHIP TO: ☐ SAME AS BILLING _____

ADDRESS: _____

SHIPPING: ☐ GROUND (FXGD) ☐ STANDARD 2 DAY (FX2D)
OVERNIGHT: ☐ PRIORITY (FX1D) ☐ 1st OVERNIGHT (FX1A)
☐ OTHER: _____

CLINICIAN: _____

CELL #: _____

PATIENT ID: _____

HEIGHT: _____ **WEIGHT:** _____ **AGE:** _____

☐ MALE ☐ FEMALE
☐ LEFT ☐ RIGHT ☐ BILATERAL

NG ENCOUNTER #: _____

MEASUREMENT DATE: _____

IN-OFFICE REQUEST DATE & TIME: _____

☐ PROJECT NEWTON (Credit applied on prior auth. denial, details on One Hanger)

HiPro is not available digitally. Please ship cast directly to HFN Houston 3620 Willowbend Blvd, Ste 1118-A | Houston, TX 77054 **HFN:** ☐ HOUSTON ☐ OTHER _____

PRE-WORK CONFIRMATION

Was the Check Brace Modified?

- ☐ No ☐ Yes (indicate modifications)*
☐ Heat Relieve ☐ Pad ☐ Trim

*

Patient was Evaluated Walking in Check Brace?

- ☐ Yes, Please Move to Definitive Device.
☐ No, Contact Me Before Moving Forward.

I am Satisfied with the Check Brace Fit (including modifications if listed).

- ☐ Yes, Please Move to Definitive Device.
☐ No, Contact Me Before Moving Forward.

BRACE DESIGN/SPECIFICATIONS

Cuff Style

- ☐ PTB ☐ Anterior Overlap
☐ Other _____

Cuff Materials

- ☐ Proflex (F9036) ☐ 1/8" Black Ice Lined

Closure Type

- ☐ Click Medical Reel ☐ Strap
☐ Other _____

Footplate

- ☐ Inner Boot ☐ 1/8" Black Ice Lined
☐ Toe Filler

Trimlines

- ☐ Malleolar Tabs ☐ Toe Lip

Finish

- ☐ Carbon ☐ Custom Fabric Provided

Click Medical Reel Location


- ☐ Posterior ☐ Lateral

REQUIREMENTS

WEIGHT

ACTIVITY LEVEL

- ☐ Low
☐ Medium
☐ High



LAB USE ONLY

Spring Category (Strut length may change due to final design)

- ☐ .5 ☐ .75 ☐ 1 ☐ 2 ☐ 3 ☐ 4 ☐ 5 ☐ 6 ☐ 7

Spring Length

- ☐ 200mm ☐ 250mm ☐ 300mm

TURNAROUND TIMES

To review current projected turnaround times for fabrication sites visit the [Daily HFN Capacity Webpage](#)