

PCC #: _____
BILL TO: _____
ADDRESS: _____

SHIP TO: SAME AS BILLING _____
ADDRESS: _____


SHIPPING: GROUND (FXGD) STANDARD 2 DAY (FX2D)
 OVERNIGHT: PRIORITY (FX1D) 1st OVERNIGHT (FX1A)
 OTHER: _____

CLINICIAN: _____
CELL #: _____
PATIENT ID: _____
HEIGHT: _____ **WEIGHT:** _____ **AGE:** _____
DIAGNOSIS: _____
AFFECTED SIDE (Check One):
 LEFT RIGHT or BILATERAL: SYMMETRICAL YES NO
ENCOUNTER #: _____
MEASUREMENT DATE: _____
IN-OFFICE REQUEST DATE & TIME: _____
 PROJECT NEWTON (Credit applied on prior auth. denial, details on One Hanger)

HFN: PHOENIX ORLANDO KANSAS HOUSTON OTHER _____


If a Discrepancy Exists, Go By Impression Measurements **Units of Measure** Millimeters Inches

PATIENT MEASUREMENTS (REQUIRED)



TUNING CRITICAL MEASUREMENTS (Optional)

Casting Block Used? Yes* No * *Best Practice: Casting block improves design accuracy and efficiency.*



Ankle Angle
 As Casted Correct to: _____°

Casting Block Setup
 Heel Height: _____ & Toe Ramp: _____

Heel Wedge
 Attached Unattached
 Shoe Heel Height = _____


Set Heel Wedge to:
 Calculate from Cast Block Setup
 Set to SVA of: _____°
 Other _____°

- **Calculated** = Casted Heel – Shoe Heel
- **SVA** = (Set AFO to SVA first) AFO Heel – Shoe Heel
- **Other** = Clinician Specified Amount

SCAN INPUT REQUIREMENTS

SCAN TYPE Split/Inside Cast Outside Cast (Preferred) **MEASUREMENTS** Average Cast Thickness _____ mm
 Positive Model: Unmodified Modified Direct Patient
 Outside Cast Forefoot ML _____ Outside Cast Ankle ML _____

.5"
.5"



Back of Heel

CLINICIAN: _____ PATIENT ID: _____

PREFERRED METHOD OF CONTACT CELL TEXT EMAIL NG MOBILE _____

DEVICE

***Weight Bearing Foot Tracing or Foot Impression Required.**

DESIGNS

Standard TMA/Toefiller

Modifications

Standard (actual foot length + 3/8") Extended Foot Length _____"

***Please provide cast in finished device alignment.
Corrections are limited to +/- 5° in the sagittal plane.**

Corrected Ankle Position

Neutral As Is Other: DF _____° PF _____°


Proximal Flare

Yes No

Sole Type (Rockers)

Ankle Joint Met-Head Heel Lisfranc MPJ Healing

(Draw key dimensions as needed)



NOTES

***Weight Bearing Foot Tracing or Foot Impression Required.
Please provide cast in finished device alignment.
Corrections are limited to +/- 5° in the sagittal plane.**

THERMOFORMING-POSTERIOR

Plastic Type

Polypropylene Copolymer

Thickness

5/32" 3/16" 1/4" Other _____

Reinforcement

Corrugation Compcore Other _____

PADDING (detail in notes section)

Alplast Plastazote Pelite/EVA PCell

Thickness: 1/8" 5/32" 3/16" 1/4"

Additional Foot Inserts: 1 2 Other _____

TRANSFER/PLASTIC COLOR Description/# _____

THERMOFORMING-ANTERIOR

Plastic Type

Polypropylene Copolymer

Thickness

5/32" 3/16" 1/4" Other _____

PADDING (detail in notes section)

Alplast Plastazote Pelite/EVA PCell

Thickness: 1/8" 5/32" 3/16" 1/4"

FINISHING

Finished Unfinished (send straps unattached)

Straps

Leave Detached Chafe Med Chafe Lat Butterfly (Double Chafe)

1" 1 1/2" 2"

Strap Material

Velcro Only Leather Back Dacron Back Other _____

Strap Color

Black White Pink Red

Beige Green Purple Blue

Non-Skid Surface

Glued: Yes No

TURNAROUND TIMES