

**PCC #:** \_\_\_\_\_

**BILL TO:** \_\_\_\_\_

**ADDRESS:** \_\_\_\_\_

**SHIP TO:**  SAME AS BILLING \_\_\_\_\_

**ADDRESS:** \_\_\_\_\_

**SHIPPING:**  GROUND (FXGD)  STANDARD 2 DAY (FX2D)  
 OVERNIGHT:  PRIORITY (FX1D)  1st OVERNIGHT (FX1A)  
 OTHER: \_\_\_\_\_

**CLINICIAN:** \_\_\_\_\_

**CELL #:** \_\_\_\_\_

**PATIENT ID:** \_\_\_\_\_

**HEIGHT:** \_\_\_\_\_ **WEIGHT:** \_\_\_\_\_ **AGE:** \_\_\_\_\_

**DIAGNOSIS:** \_\_\_\_\_

**AFFECTED SIDE (Check One)**

LEFT  RIGHT or  BILATERAL: SYMMETRICAL  YES  NO

**ENCOUNTER #:** \_\_\_\_\_


**MEASUREMENT DATE:** \_\_\_\_\_

**IN-OFFICE REQUEST DATE & TIME:** \_\_\_\_\_

**HFN:**  KANSAS  ORLANDO  PHOENIX


**If a Discrepancy Exists, Go By:**  Impression  Measurements **Units of Measure:**  Millimeters  Inches

**PATIENT MEASUREMENTS (REQUIRED)**



**ALIGNMENT CASTING BLOCK/TUNING (Optional)**

**Alignment Casting Block Used?**  Yes\*  No \*Best Practice: Casting block improves design accuracy and efficiency.



**Ankle Angle**

As Casted  Correct to: \_\_\_\_\_°

**Casted on**

Cast Block Height: \_\_\_\_\_ & Toe Ramp: \_\_\_\_\_

**External Heel Wedge**

Attached  Unattached

Shoe Heel Height = \_\_\_\_\_

**Set Heel Wedge to**

Calculate from Cast Block Setup

Set to SVA of: \_\_\_\_\_°

Other: \_\_\_\_\_°

• **Calculated** = Casted Heel – Shoe Heel

• **SVA** = (Set AFO to SVA first) AFO Heel – Shoe Heel

• **Other** = Clinician Specified Amount

**DIGITAL SCAN INPUT REQUIREMENTS**

**SCAN TYPE:**  Split/Inside Cast  Outside Cast (Preferred)  
 Positive Model:  Unmodified  Modified  Direct Patient

**MEASUREMENTS:** Average Cast Thickness \_\_\_\_\_ mm  
 Outside Cast Forefoot ML \_\_\_\_\_ Outside Cast Ankle ML \_\_\_\_\_

CLINICIAN: \_\_\_\_\_ PATIENT ID: \_\_\_\_\_

PREFERRED METHOD OF CONTACT:  CELL  TEXT  EMAIL  MICROSOFT TEAMS \_\_\_\_\_

**DEVICE**

**Options**

Single Upright  Double Upright (CORRECTED TRACING IS REQUIRED)

**Modifications**

Standard  Tone Reduction  ST Mod

Intrinsic Heel Mod \_\_\_\_\_°  Medial  Lateral

Additional Build Ups/Reductions

Heel Height:  None  Other \_\_\_\_\_

**Corrected Ankle Position**

Neutral  As Is Other:  DF \_\_\_\_\_°  PF \_\_\_\_\_°

**Final Corrected Forefoot Position**

Right:  Neutral  As Is  Other \_\_\_\_\_

Left:  Neutral  As Is  Other \_\_\_\_\_

**Final Corrected Hindfoot Position**

Right:  Neutral  As Is  Other \_\_\_\_\_

Left:  Neutral  As Is  Other \_\_\_\_\_

**DESIGN**

**Material**

AL  SS  TI

**Bar Size**

1/4 x 3/4"  1/4 x 5/8"  3/16 x 1/2"  3/16 x 3/4"

3/16 x 5/8"  1/8 x 1/2"  Other \_\_\_\_\_

**Pre-Set**

Pins  Springs  DF \_\_\_\_\_°  PF \_\_\_\_\_°

**Stirrup**

Solid  Solid-Extended  Solid-White  Split  UCBL

Other \_\_\_\_\_

**Ankle Joints**

DAAJ  Free Motion  Limited Action  Dorsi Assist

Other \_\_\_\_\_

**DESIGN (continued)**

**Knee Joints**

Bail Lock  Drop Lock  Cam Lock  Dial Lock

Step Lock  Spring Lever Lock  Ball Retainers

Trigger Release  Posterior Offset  Polycentric

**Finish Options**

Growth Extensions  Tibial Torsion Spreader Bar \_\_\_\_\_°

**FINISHING**

**Straps**

Eyelets  Calf Lacer  Thigh Lacer  Instep

PTB  POS Check  4-Buckle  5-Buckle

T Strap  Medial  Lateral

Other \_\_\_\_\_

**Strap Material**

Leather Back (F0046)  Dacron Back (F0046)  Other \_\_\_\_\_

**Strap Color**

Black  White  Pink  Red

Beige  Green  Purple  Blue

**Shoes**

Heel Lift \_\_\_\_\_

Heel & Sole Lift \_\_\_\_\_

Medial Wedge \_\_\_\_\_

Lateral Wedge \_\_\_\_\_

Steel Shank \_\_\_\_\_

Other \_\_\_\_\_

**NOTES**

Large empty box for notes.

**TURNAROUND TIMES**

To review current projected turnaround times for fabrication sites visit the [Daily HFN Capacity Webpage](#).