

Pediatric



ORTHOTICS Master the Catwalk

Meerkat SCOOTZ

PEDIATRIC LOW TONE SUPRAMALLEOLAR ORTHOSIS (SMO)

Developed by Hanger Clinic clinicians, engineers, and technicians, Meerkats are designed to accelerate each child's specific rehabilitation goals. Individual patient measurements and data-controlled inputs are used to guarantee fit. Our commitment is to provide quality, best value products so practitioners can concentrate on doing what they do best, focused patient care.

INDICATIONS

- Hypotonia
- Flexible pes planus – moderate to severe
- Ankle instability – moderate to severe

CONTRAINDICATIONS

- Clubfoot
- Forefoot deformities
- High tone

ENHANCED FEATURES

- Sustentaculum tali modification
- Medial external and intrinsic heel posting
- No cast required for children <80 lbs



Did you know?

- Meerkats live in South Africa.
- They live in groups called “mobs” which consist of 2-3 families.
- Standing on their hind legs helps them survey the environment and absorb heat on their bellies after a long, cold desert night.

