

# HANGER TMAs

Transmetatarsal amputations

These TMAs were developed by a team of highly skilled Hanger Clinic orthetists who combined their clinical experience, knowledge of anatomy, and computerized gait analysis to completely rethink and improve designs.



## TMA-1 Indications

- Fixed/Painful Ankle

## Design

- Mid-Calf
- Copoly
- Solid AFO
- Northvane SMO
- TMA Insert



## TMA-2 Indications

- Compromised Skin
- Equino-Varus
- PF Contracture

## Design

- Mid-Calf
- Copoly
- Articulated AFO
- Northvane SMO
- TMA Insert
- Launchpad Joint
- X