

This is not the work order. **DO NOT SUBMIT THIS SHEET.** After recording your measurements complete the online work order form at hangerfabrication.com or scan the QR code to the right.



PATIENT ID: (first initial) _____ (last name) _____

MEASUREMENTS

***IMPORTANT:** Provide trim length from **Distal Limb**.
Fabrication will adjust length to accommodate for end pad.

The diagram illustrates the measurement process for a custom transfemoral amputee. It features a central shaded area representing the limb, with various measurement points and labels:

- Left and Right Limb Markers:** Two columns of circles on the left and right sides, labeled "Left" and "Right" respectively. The circles are numbered from "0 Lvl" down to "-12" in increments of 2 inches.
- Zero Level at Medial Trim Height:** A horizontal line with an arrow pointing to the "0 Lvl" circle, indicating the starting point for measurements.
- Medial Limb Trim Height:** Two boxes labeled "Left" and "Right" for recording measurements at the medial limb trim height.
- Lateral Limb Trim Height:** Two boxes labeled "Left" and "Right" for recording measurements at the lateral limb trim height.
- Distal Limb:** A vertical dashed line extending from the "Zero Level" down to the bottom of the limb, indicating the distal limb measurement point.

NOTES