

LOW-TONE SMO

Supramalleolar Orthosis

Developed by Hanger Clinic clinicians, engineers, and technicians, these SMOs are designed to accelerate each child's specific rehabilitation goals. Individual patient measurements and data-controlled inputs are used to guarantee fit. Our commitment is to provide quality, best value products so practitioners can concentrate on doing what they do best, focused patient care.



Features

- Sustentaculum tali modification
- Medial external and intrinsic heel posting
- No cast required for children < 80 lbs.

Indications

- Hypotonia
- Flexible pes planus – moderate to severe
- Ankle instability – moderate to severe

Contraindications

- Clubfoot
- Forefoot deformities
- High tone

