




# Partial Hand

## Frequently Asked Questions

### What's the process for Partial Hand Delivery?



Please refer to **Partial Hand Clinician Steps** for more details.

### How much does a Partial Hand cost?

#### Can I get a quote?

We are pleased to provide quotes on request. Each partial hand is designed uniquely for each patient. Note: we can provide estimations on device FABRICATION only, componentry costs are not included. In order to provide a quote, please provide us with amputation level and number of fingers to be fabricated.

#### Do I need to order parts?

Yes. Clinicians are responsible for ordering components and providing them along with the device design specifications.

#### Where can I order the Partial Hand parts?

**Hanger Direct, (internal) – Partial Hand and Finger 3D Print Catalog 2022** or College Park (external)

#### How do I cast the patient?

We recommend using Body Double or P386 (Gray) impression material. Cast requirements are dependent on the type of device you are ordering:

- Full M-Fingers: up to the cubital fold
- Partial M-Fingers: 6" past the wrist
- Titans: 2" past wrist
- Cast should be 1/2" thick and send photos as well (*Best practice: Sound side mold is not required but ok to send as well.*)

#### How long does it take to finish?

Device fabrication time can vary with the design. In general:

- Silicone only requires less time (2-3 weeks).  
Initial Fitting 2-3 weeks
- Fitting devices with a prepreg frame takes approximately 3-5 weeks.
- Final device: 2-4 weeks from date of return.  
Final Fitting 2-4 weeks.

Fitting dates cannot be determined until all parts to support the project have been received. Parts can take up to 2 weeks to arrive.

#### What determines the device features?

Each device is designed to promote patient functionality. Please contact us to discuss design needs and current options that are available.